

Laboratory Test Updates

(effective on January 09, 2023)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	Update In Additional	New Tests	Internalized Tests	Reactivated Tests	Inactive
1	C8301839	s-Phenylmercapturic Acid							X			
2	C8208501	Aldolase								X		
3	C8557602	Platelet Function Test		X	X	X		X				
4	C8259501	Cryoglobulin				X						
5	C8524002	ADAM TS 13 (von Willebrand Factor Splitting Protease)			X							

1	s-Phenylmercapturic Acid Test Code : 8301839 - Laboratory Services			
	Methodology : LC-MS/MS - Performed : Mon-Fri 08:00 - Reported : 10 days			
	Specimen Type	Container	Volume	Transport
	Urine	Plastic Collection Cup	3mL	2-8°C

2	Aldolase Test Code : 8208501 - Laboratory Services			
	Methodology : ENZYMATIC - Performed : Tue 08:00 - Reported : 2 weeks			
	Specimen Type	Container	Volume	Transport
	Serum	Separator Tube	1mL	2-8°C

3	Platelet Function Test Test Code : 8557602 - Laboratory Services			
	Methodology : IMPEDANCE AGGREGOMETRY - Performed : Mon-Fri 16:00 Sat 12:00 - Reported : Tue-Sat 18:00 Mon 18:00			
	Specimen Type	Container	Volume	Transport
	Whole Blood	Hirudin Tube	2mL	18-24°C
	Additional Information ▼			
	Whole Blood Programmed testing! For the platelet function test, the blood is collected into anticoagulant tubes with hirudin. Specimen should be sent to the laboratory immediately, to ensure testing within 3 hours of collection. For further information, contact the laboratory.			

4	Cryoglobulin Test Code : 8259501 - Laboratory Services			
	Methodology : MANUAL - Performed : Mon-Fri - Reported : 1 weeks			
	Specimen Type	Container	Volume	Transport
	Special sample type, Apply to Laboratory	Plain Tube		37°C
Additional Information ▼				
<p>Special sample type, Apply to Laboratory</p> <p>The blood specimen must be collected in two 10 mL warm tubes (37°C) in the absence of anticoagulants. The sample should be immediately transferred to a 37°C incubator/water-bath and allowed to clot and settle at 37°C for an for 1 hr (DO NOT centrifuge). The self-separated serums are transferred to 2 secondary tubes of at least 5 mL in each tube and transported to the laboratory at 37°C.</p>				

5	ADAM TS 13 (von Willebrand Factor Splitting Protease) Test Code : 8524002 - Laboratory Services			
	Methodology : ELISA - Performed : Mon 08:00 - Reported : 2 weeks 18:00			
	Specimen Type	Container	Volume	Transport
Plasma (Citrate), Frozen	Blue (Sodium citrate)	4mL	<0°C	