

Adrenal Stress Test

Test Code : **8460005** - Laboratory Services

1

Methodology : **ECLIA, LC-MS/MS** - Performed : **Wed 08:00** - Reported : **Thu 18:00**

Specimen Type	Container	Volume	Transport
Apply to Laboratory	Call the Laboratory		<0°C

Additional Information

Apply to Laboratory

4 Saliva samples should be collected and their hours are as follows : Saliva morning (07:00-09:00), Saliva noon (11:00-13:00), Saliva afternoon (16:00-18:00) and Saliva midnight (22:00-00:00). Considerations to be taken before the sample is collected : Do not brush teeth before collecting specimen. Do not eat or drink for 15 minutes prior to specimen collection. Record the exact collection date and time on both the container label on the swab tube. To Use the Salivette : Remove top cap of the tube to expose swap from the salivette. Place swap directly into the mouth by tipping the tube so the swap falls into the mouth. Keep the swap in the mouth and chew it gently for about 2 minutes. Spit the swap back into the tube. Replace the cap. Make sure cap is pushed on tightly. Record the exact collection date and time on the tube. Along with the samples, the Adrenal Stress Test Information Form should also be filled out and sent.

Neuron Specific Enolase (NSE)

Test Code : **8631605** - Laboratory Services

3

Methodology : **ECLIA** - Performed : **Mon 08:00 Fri 08:00** - Reported : **Mon 18:00 Mon 18:00**

Specimen Type	Container	Volume	Transport
Serum	Separator Tube	2mL	2-8°C

Neuron Specific Enolase (NSE), CSF

Test Code : **8631610** - Laboratory Services

4

Methodology : **ECLIA** - Performed : **Mon 08:00 Fri 08:00** - Reported : **Mon 18:00 Mon 18:00**

Specimen Type	Container	Volume	Transport
CSF	Sterile CSF Tube	1mL	2-8°C

CA 72-4Test Code : **8631604** - Laboratory Services**5**Methodology : **ECLIA** - Performed : **Mon 08:00 Fri 08:00** - Reported : **Mon 18:00 Mon 18:00**

Specimen Type	Container	Volume	Transport
Serum	Separator Tube	1mL	2-8°C

Complement C1qTest Code : **8616011** - Laboratory Services**6**Methodology : **EIA** - Performed : **Mon-Thu 08:00 Fri 08:00** - Reported : **Tue-Fri 18:00 Mon 18:00**

Specimen Type	Container	Volume	Transport
Serum	Separator Tube	1mL	2-8°C