

Laboratory Test Updates

(effective on March 11, 2024)

| | Test Code | Test Name | Name Change | Methodology | Performed/Reported Schedule | Specimen Requirements | Reference Interval | New Tests | Internalized Tests | Reactivated Tests | Inactive |
|---|-----------|---|-------------|-------------|-----------------------------|-----------------------|--------------------|-----------|--------------------|-------------------|----------|
| 1 | C8414300 | 17-hydroxypregnenolon, Urine (24 hours) | | | | | | | | | X |
| 2 | C8294301 | Glucagon | | | | X | | | | | |

| | | | | |
|---|---|----------------------------|--------|-----------|
| 2 | <p>Glucagon</p> <p>Test Code : 8294301 - Laboratory Services</p> | | | |
| | Methodology : RIA - Performed : Tue 08:00 - Reported : 17 days | | | |
| | Specimen Type | Container | Volume | Transport |
| | Plasma | Lavender (EDTA) (Trasylol) | 2mL | <0°C |
| | Additional Information ▼ | | | |
| <p>Plasma</p> <p>Patient should fast overnight. Obtain EDTA-Trasylol (APROTININ) tube from laboratory. Sample should be centrifuged within 30 mins of collection. After centrifugation, the plasma is separated and transferred to the laboratory.</p> | | | | |

