

Laboratory Test Updates

(effective on May 12, 2025)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8818077	Daratumumab Interference Test (MALDI-TOF)										X
2	C8538501	Fibrinogen				X						
3	C8301809	Short Chain Fatty Acids, Stool						X				
4	C8019708	Specific IgE Inhalant Panel (26 parameters)			X							
5	C8019709	Specific IgE Food Allergen (34 parameters)			X							
6	C8378907	Silver	X			X	X					

2	Fibrinogen Test Code : 8538501 - Laboratory Services		
	Methodology : COAGULOMETRY - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
	Plasma (Citrate)	Blue (Sodium citrate)	2mL
	Additional Information		▼
	The plasma sample should be frozen by portioning it into at least two separate tubes.		

3	Short Chain Fatty Acids, Stool Test Code : 8301809 - Laboratory Services		
	Method : GC - Performed : Fri 08:00 - Reported : 15 days		
	Sample Name	Collection Container/Tube	Specimen Volume
	Stool, Frozen	Plastic Collection Cup	10gr
	Transfer		<0°C

Additional Information

This panel includes: Acetic Acid, Propionic Acid, I-Butyric Acid, N-Butyric Acid, I-Valeric Acid, N-Valeric Acid, N-Caproic Acid ve N-Heptanoic Acid.

4	Specific IgE Inhalant Panel (26 parameters) Test Code : 8019708 - Laboratory Services		
	Method : IMMUNOBLOTTING - Performed : Mon 08:00 Thu 08:00 - Reported : Wed 18:00 Mon 18:00		
	Sample Name	Collection Container/Tube	Specimen Volume
	Serum	Separator Tube	0,5mL
	Transfer		2-8°C

5	Specific IgE Food Allergen (34 parameters) Test Code : 8019709 - Laboratory Services		
	Methodology : IMMUNOBLOTTING - Performed : Mon 08:00 Thu 08:00 - Reported : Wed 18:00 Mon 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
	Serum	Separator Tube	0,5mL
	Transport		2-8°C

6	Silver Test Code : 8378907 - Laboratory Services		
	Method : ICP-MS - Performed : Tue 08:00 - Reported : 2 weeks		
	Sample Name	Collection Container/Tube	Specimen Volume
	Plasma (Heparin)	Sodium Heparin-Royal Blue	2mL
	Transfer		2-8°C

Reference Ranges

Normal Values		
Gender	Lower Value	Unit
	<0,3	ug/L