

Laboratory Test Updates

(effective on Jan 17, 2022)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	Update In Additional	New Tests	Internalized Tests	Reactivated Tests	Inactive
1	C8349806	Troponin I-High Sensitive						X				
2	C8388301	Cystatine-C			X							
3	C8388303	Cysteine-C, Urine			X							
4	C8616011	Complement C1q					X					

1	Troponin I-High Sensitive			
	Test Code : 8349806 - Laboratory Services			
	Method : CLIA - Performed : Mon-Sat - Reported : 2 hours			
	Sample Name	Container	Amount	Transfer
Serum	Separator Tube	4mL	2-8°C	
Plasma (Heparin)	Heparin-Yeşil Kapak(Jelli)	4mL	2-8°C	

2	Cystatine-C			
	Test Code : 8388301 - Laboratory Services			
	Method : NEPHELOMETRY - Performed : Mon 08:00 - Reported : 15 days			
Sample Name	Container	Amount	Transfer	
Serum, Frozen	Separator Tube	1,5mL	<0°C	

3	Cysteine-C, Urine			
	Test Code : 8388303 - Laboratory Services			
	Method : NEPHELOMETRY - Performed : Mon 08:00 - Reported : 15 days			
Sample Name	Container	Amount	Transfer	
Urine, Frozen	Plastic Collection Cup	2mL	<0°C	

4	Complement C1q			
	Test Code : 8616011 - Laboratory Services			
	Method : NEPHELOMETRY - Performed : Tue 08:00 - Reported : 1 weeks			
	Sample Name	Container	Amount	Transfer
	Serum	Separator Tube	1mL	2-8°C
	Reference Ranges ▼			
Normal Values				
Gender	Lower Value	Upper Value	Unit	
	15.7	30.6	mg/dL	

--	--