



| 1             | C-Peptide  |           |        |           |
|---------------|--|-----------|--------|-----------|
|               | Test Code : <b>8468101</b> - Laboratory Services   |           |        |           |
|               | Methodology : <b>CLIA</b> - Performed : <b>Mon-Sat 12:00</b> - Reported : <b>Mon-Sat 18:00</b> |           |        |           |
|               | Specimen Type  | Container | Volume | Transport |
| Serum, Frozen | Separator Tube   | 1mL       | <0°C   |           |

| 2 | C-Peptide, Postprandial  |           |        |           |
|---|--|-----------|--------|-----------|
|   | Test Code : <b>8468103</b> - Laboratory Services   |           |        |           |
|   | Methodology : <b>CLIA</b> - Performed : <b>Mon-Sat 12:00</b> - Reported : <b>Mon-Sat 18:00</b> |           |        |           |
|   | Specimen Type  | Container | Volume | Transport |
|   |  |           |        |           |

| 3     | Cryptococcus neoformans antigen, Serum   |           |        |           |
|-------|--|-----------|--------|-----------|
|       | Test Code : <b>8732699</b> - Laboratory Services   |           |        |           |
|       | Methodology : <b>IMMUNOCHROMATOGRAPHY</b> - Performed : <b>Mon-Sat 08:00</b> - Reported : <b>Mon-Sat 18:00</b> |           |        |           |
|       | Specimen Type  | Container | Volume | Transport |
| Serum | Separator Tube   | 1mL       | 2-8°C  |           |

| 4   | Cryptococcus neoformans antigen, BOS   |           |        |           |
|-----|--|-----------|--------|-----------|
|     | Test Code : <b>8732700</b> - Laboratory Services   |           |        |           |
|     | Methodology : <b>IMMUNOCHROMATOGRAPHY</b> - Performed : <b>Mon-Sat 08:00</b> - Reported : <b>Mon-Sat 18:00</b> |           |        |           |
|     | Specimen Type  | Container | Volume | Transport |
| CSF | Biologic Fluid Tube  | 1mL       | 2-8°C  |           |