



Laboratory Test Updates

(effective on September 01, 2022)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	Update In Additional	New Tests	Internalized Tests	Reactivated Tests	Inactive
1	C8616101	C1 Esterase Inhibitor, Antigen					X					
2	C8666502	Epstein-Barr Virus (EBV)-VCA IgG		X								
3	C8267901	Estron (E1)		X			X			X		
4	C8625619	Thyroid Stimulan Immunglobulin (TSI)		X								
5	C8223201	Beta-2-Microglobulin, Serum			X							
6	C8223203	Beta-2-Microglobulin,urine(spot)			X							
7	C8223204	Beta-2-Microglobulin,urine(24 hr)			X							

1	C1 Esterase Inhibitor, Antigen			
	Test Code : 8616101 - Laboratory Services			
	Method : TURBIDIMETRY - Performed : Thu 08:00 - Reported : Tue 18:00			
	Sample Name	Container	Amount	Transfer
	Serum	Separator Tube	2mL	2-8°C
	Reference Ranges 			
Normal Values				
Gender	Lower Value	Upper Value	Unit	
	21	38	mg/dL	

2	Epstein-Barr Virus (EBV)-VCA IgG			
	Test Code : 8666502 - Laboratory Services			
	Method : CMIA - Performed : Mon-Thu 08:00 Fri 08:00 - Reported : Tue-Fri 18:00 Mon 18:00			
	Sample Name	Container	Amount	Transfer
Serum	Separator Tube	1mL	2-8°C	

3	Estron (E1)				
	Test Code : 8267901 - Laboratory Services				
	Method : ELISA - Performed : Tue 08:00 - Reported : 17 days				
	Sample Name	Container	Amount	Transfer	
	Serum	Separator Tube	0,3mL	2-8°C	
	Reference Ranges 				
	Normal Values				
	Gender	Age Range	Lower Value	Upper Value	Unit
	Male	13-100 YEAR	0	174	pg/mL
	Special Values				
Gender	Description	Age Range	Lower Value	Upper Value	Unit
Female	Follicular Phase	12-100 YEAR	18.9	176.1	pg/mL
Female	Luteal Phase	12-100 YEAR	37.1	199.8	pg/mL
Female	Postmenopause	12-100 YEAR	0	125.1	pg/mL

4	Thyroid Stimulan Immunglobulin (TSI)			
	Test Code : 8625619 - Laboratory Services			
	Method : CLIA - Performed : Mon-Thu 08:00 Fri 08:00 - Reported : Tue-Fri 18:00 Mon 18:00			
	Sample Name	Container	Amount	Transfer
Serum	Separator Tube	0,5mL	2-8°C	

5	Beta-2-Microglobulin, Serum			
	Test Code : 8223201 - Laboratory Services			
	Method : NEPHELOMETRY - Performed : Fri 08:00 - Reported : Mon 18:00			
	Sample Name	Container	Amount	Transfer
Serum	Separator Tube	0,5mL	2-8°C	

6	Beta-2-Microglobulin,urine(spot)			
	Test Code : 8223203 - Laboratory Services			
	Method : NEPHELOMETRY - Performed : Fri 08:00 - Reported : Mon 18:00			
	Sample Name	Container	Amount	Transfer
Urine, Random	Plastic Collection Cup	5mL	2-8°C	

7	Beta-2-Microglobulin,urine(24 hr)			
	Test Code : 8223204 - Laboratory Services			
	Method : NEPHELOMETRY - Performed : Fri 08:00 - Reported : Mon 18:00			
	Sample Name	Container	Amount	Transfer
24-hour Urine	Plastic 24-Hour Urine Container	10mL	2-8°C	