

Laboratory Test Updates

(effective on August 21, 2023)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	Update In Additional	New Tests	Internalized Tests	Reactivated Tests	Inactive
1	C8255401	Creatin Kinase Isoenzymes						X				
2	C8371501	Lipoprotein Electrophoresis			X							
3	C8440302	Testosterone, Bioavailable			X							
4	C8294101	Gastrin			X							
5	C8265708	Alpha-Galactosidase, Leucocyte			X							
6	C8265698	Alpha-Fucosidase, Leukocytes			X							
7	C8265711	Alpha-Glucosidase, Leucocyte			X							
8	C8265697	Alpha-Iduronidase, Leukocytes			X							
9	C8265699	Alpha-Mannosidase, Leukocytes			X							
10	C8265704	Arylsulphatase A, Leucocyte			X							
11	C8265703	Beta-galactosidase, Leucocyte			X							
12	C8296302	Beta-glucosidase, Leucocyte			X							
13	C8308001	Beta-Hexosaminidase A			X							
14	C8308002	Beta-Hexosaminidase B			X							
15	C8265706	Galactocerebrosidase, Leucocyte			X							
16	C8265705	Sphingomyelinase, Leucocyte			X							

17	C8416502	Protein Electrophoresis, Urine			X								
18	C8408001	Alkaline Phosphatase, Isoenzyme Electrophoresis			X								
19	C8278405	IgA Subgroups			X								
20	C8278404	IgD			X								
21	C8362501	Lactate Dehydrogenase (LDH) Isoenzymes			X								
22	C8388301	Cystatine-C			X								
23	C8388303	Cysteine-C, Urine			X								
24	C8744901	Aspergillus Galactomannan Antigen			X	X							

1	Creatin Kinase Isoenzymes Test Code : 8255401 - Laboratory Services			
	Methodology : ELECTROPHORESIS - Performed : Mon-Fri 08:00 - Reported : 1 week			
	Specimen Type	Container	Volume	Transport
	Serum	Separator Tube	1mL	2-8°C
	Additional Information ▼			
Hemolysed sample is not acceptable.				

2	Lipoprotein Electrophoresis Test Code : 8371501 - Laboratory Services			
	Methodology : ELECTROPHORESIS - Performed : Wed 08:00 - Reported : Fri 18:00			
	Specimen Type	Container	Volume	Transport
	Serum	Separator Tube	1mL	2-8°C

3	Testosterone, Bioavailable Test Code : 8440302 - Laboratory Services			
	Methodology : - Performed : Mon 08:00 Thu 08:00 - Reported : Mon 18:00 Thu 18:00			
	Specimen Type	Container	Volume	Transport
	Serum	Separator Tube	1mL	2-8°C

4	Gastrin Test Code : 8294101 - Laboratory Services			
	Methodology : CLIA - Performed : Wed 08:00 - Reported : Thu 18:00			
	Specimen Type	Container	Volume	Transport
	Serum, Frozen	Separator Tube	1mL	<0°C

5	Alpha-Galactosidase, Leucocyte Test Code : 8265708 - Laboratory Services			
	Methodology : FLUOROMETRY - Performed : Thu 08:00 - Reported : Tue 18:00			
	Specimen Type	Container	Volume	Transport
	Whole Blood (EDTA)	Lavender (EDTA)	8mL	2-8°C
	Whole Blood (ACD)	ACD-Yellow Top	8mL	2-8°C

6	Alpha-Fucosidase, Leukocytes Test Code : 8265698 - Laboratory Services			
	Methodology : FLUOROMETRY - Performed : Thu 08:00 - Reported : Tue 18:00			
	Specimen Type	Container	Volume	Transport
	Whole Blood (EDTA)	Lavender (EDTA)	20mL	2-8°C
	Whole Blood (ACD)	ACD-Yellow Top	8mL	2-8°C

7	Alpha-Glucosidase, Leucocyte Test Code : 8265711 - Laboratory Services			
	Methodology : ENZYMATIC - Performed : Thu 08:00 - Reported : Tue 18:00			
	Specimen Type	Container	Volume	Transport
	Whole Blood (EDTA)	Lavender (EDTA)	8mL	2-8°C
	Whole Blood (ACD)	ACD-Yellow Top	8mL	2-8°C

8	Alpha-Iduronidase, Leukocytes Test Code : 8265697 - Laboratory Services			
	Methodology : FLUOROMETRY - Performed : Thu 08:00 - Reported : Tue 18:00			
	Specimen Type	Container	Volume	Transport
	Whole Blood (EDTA)	Lavender (EDTA)	20mL	2-8°C
	Whole Blood (ACD)	ACD-Yellow Top	8mL	2-8°C

9	Alpha-Mannosidase, Leukocytes Test Code : 8265699 - Laboratory Services			
	Methodology : FLUOROMETRY - Performed : Thu 08:00 - Reported : Tue 18:00			
	Specimen Type	Container	Volume	Transport
	Whole Blood (EDTA)	Lavender (EDTA)	20mL	2-8°C
Whole Blood (ACD)	ACD-Yellow Top	8mL	2-8°C	

10	Arylsulphatase A, Leucocyte Test Code : 8265704 - Laboratory Services			
	Methodology : PHOTOMETRY - Performed : Thu 08:00 - Reported : Tue 18:00			
	Specimen Type	Container	Volume	Transport
	Whole Blood (EDTA)	Lavender (EDTA)	10mL	2-8°C
Whole Blood (ACD)	ACD-Yellow Top	8mL	2-8°C	

11	Beta-galactosidase, Leucocyte Test Code : 8265703 - Laboratory Services			
	Methodology : FLUOROMETRY - Performed : Thu 08:00 - Reported : Tue 18:00			
	Specimen Type	Container	Volume	Transport
	Whole Blood (EDTA)	Lavender (EDTA)	20mL	2-8°C
Whole Blood (ACD)	ACD-Yellow Top	8mL	2-8°C	

12	Beta-glucosidase, Leucocyte Test Code : 8296302 - Laboratory Services			
	Methodology : FLUOROMETRY - Performed : Thu 08:00 - Reported : Tue 18:00			
	Specimen Type	Container	Volume	Transport
	Whole Blood (EDTA)	Lavender (EDTA)	8mL	2-8°C
Whole Blood (ACD)	ACD-Yellow Top	8mL	2-8°C	

13	Beta-Hexosaminidase A Test Code : 8308001 - Laboratory Services			
	Methodology : FLUOROMETRY - Performed : Thu 08:00 - Reported : Tue 18:00			
	Specimen Type	Container	Volume	Transport
	Plasma (EDTA)	Lavender (EDTA)	1mL	2-8°C
	Plasma	ACD-Yellow Top	1mL	2-8°C

14	Beta-Hexosaminidase B Test Code : 8308002 - Laboratory Services			
	Methodology : FLUOROMETRY - Performed : Thu 08:00 - Reported : Tue 18:00			
	Specimen Type	Container	Volume	Transport
	Whole Blood (EDTA)	Lavender (EDTA)	10mL	2-8°C
	Whole Blood (ACD)	ACD-Yellow Top	8mL	2-8°C

15	Galactocerebrosidase, Leucocyte Test Code : 8265706 - Laboratory Services			
	Methodology : FLUOROMETRY - Performed : Thu 08:00 - Reported : Tue 18:00			
	Specimen Type	Container	Volume	Transport
	Whole Blood (EDTA)	Lavender (EDTA)	8mL	2-8°C
	Whole Blood (ACD)	ACD-Yellow Top	8mL	2-8°C

16	Sphingomyelinase, Leucocyte Test Code : 8265705 - Laboratory Services			
	Methodology : FLUOROMETRY - Performed : Thu 08:00 - Reported : Tue 18:00			
	Specimen Type	Container	Volume	Transport
	Whole Blood (EDTA)	Lavender (EDTA)	8mL	2-8°C
	Whole Blood (ACD)	ACD-Yellow Top	8mL	2-8°C

17	Protein Electrophoresis, Urine Test Code : 8416502 - Laboratory Services			
	Methodology : ELECTROPHORESIS - Performed : Tue 08:00 Fri 08:00 - Reported : Fri 18:00 Wed 18:00			
	Specimen Type	Container	Volume	Transport
24-hour Urine	Plastic 24-Hour Urine Container	10mL	2-8°C	

18	Alkaline Phosphatase, Isoenzyme Electrophoresis Test Code : 8408001 - Laboratory Services			
	Methodology : ELECTROPHORESIS - Performed : Wed 08:00 - Reported : 1 week			
	Specimen Type	Container	Volume	Transport
Serum	Separator Tube	1mL	2-8°C	

19	IgA Subgroups Test Code : 8278405 - Laboratory Services			
	Methodology : TURBIDIMETRY - Performed : Thu 08:00 - Reported : 15 days			
	Specimen Type	Container	Volume	Transport
Serum, Frozen	Separator Tube	1mL	<0°C	

20	IgD Test Code : 8278404 - Laboratory Services			
	Methodology : TURBIDIMETRY - Performed : Thu 08:00 - Reported : Fri 18:00			
	Specimen Type	Container	Volume	Transport
Serum, Frozen	Separator Tube	1mL	<0°C	

21	Lactate Dehydrogenase (LDH) Isoenzymes Test Code : 8362501 - Laboratory Services			
	Methodology : ELECTROPHORESIS - Performed : Tue 08:00 - Reported : 1 week			
	Specimen Type	Container	Volume	Transport
Serum	Separator Tube	1mL	2-8°C	

22	Cystatine-C Test Code : 8388301 - Laboratory Services			
	Methodology : NEPHELOMETRY - Performed : Tue 08:00 - Reported : Wed 18:00			
	Specimen Type	Container	Volume	Transport
	Serum, Frozen	Separator Tube	1,5mL	<0°C

23	Cysteine-C, Urine Test Code : 8388303 - Laboratory Services			
	Methodology : NEPHELOMETRY - Performed : Tue 08:00 - Reported : Wed 18:00			
	Specimen Type	Container	Volume	Transport
	Urine, Frozen	Plastic Collection Cup	2mL	<0°C

24	Aspergillus Galactomannan Antigen Test Code : 8744901 - Laboratory Services			
	Methodology : CLIA - Performed : Mon-Thu 08:00 Fri 08:00 - Reported : Tue-Fri 18:00 Mon 18:00			
	Specimen Type	Container	Volume	Transport
	Bronchoalveolar Lavage	Biologic Fluid Tube	1mL	2-8°C
	Serum	Separator Tube	1mL	2-8°C
	Plasma (EDTA)	Lavender (EDTA)	1mL	2-8°C
	CSF	Sterile CSF Tube	1mL	2-8°C