



1	<b>NMO-IgG Antibody (Anti-Aquaporin 4 Antibody), CSF (IFA)</b> Test Code : <b>8625626</b> - Laboratory Services			
	Methodology : <b>IFT</b> - Performed : <b>Tue 08:00</b> - Reported : <b>Wed 18:00</b>			
	Specimen Type	Container	Volume	Transport
CSF	Sterile CSF Tube	1mL	2-8°C	

2	<b>NMO-IgG Antibody (Anti-Aquaporin 4 Antibody) (IFA)</b> Test Code : <b>8625620</b> - Laboratory Services			
	Methodology : <b>IFT</b> - Performed : <b>Tue 08:00</b> - Reported : <b>Wed 18:00</b>			
	Specimen Type	Container	Volume	Transport
Serum	Separator Tube	2mL	2-8°C	

4	<b>Quad Combo Test (RSV-Adenovirus Antigen /Influenza (A+B) Antigen/ SARS-CoV-2 Ag)</b> Test Code : <b>8392916</b> - Laboratory Services			
	Methodology : <b>IMMUNOCHROMATOGRAPHY</b> - Performed : <b>Mon-Sun</b> - Reported : <b>1 hours</b>			
	Specimen Type	Container	Volume	Transport
Nasopharyngeal Swab	Call the Laboratory		2-8°C	

## Clonazepam


Test Code : **8029916** - Laboratory Services

5

Methodology : **HPLC** - Performed : **Tue 08:00** - Reported : **1 weeks**

Specimen Type	Container	Volume	Transport
Serum	Plain Tube	1mL	2-8°C

Additional Information 

Reference Ranges 

### Special Values

Gender	Description	Lower Value	Upper Value	Unit
	Anticonvulsant	20	70	ng/mL
	Sleep disorder/Anxiolytic	4	80	ng/mL
	Toksik düzey		>80	ng/mL