

Modified acid fast stain

Test Code : **8720602** - Laboratory Services

1

Method : **MANUAL** - Performed : **Mon-Sat 08:00** - Reported : **Mon-Sat 18:00**

Sample Name	Container	Amount	Transfer
Stool	Plastic Collection Cup	5gr	2-8°C
Sputum	Plastic Collection Cup		2-8°C
Bronchoalveolar Lavage	Biologic Fluid Tube		2-8°C
Transtracheal Aspirate	Biologic Fluid Tube		2-8°C
Pus	Biologic Fluid Tube		2-8°C
Other	Sterile Tube		2-8°C

Additional Information

Stool should be collected prior to the use of antibiotics or antidiarrheal medications. Sample should be given at least 48 hours after endoscopy, colonoscopy procedures. Sample should be collected without contact with urine or toilet water.

Mevalonic Acid, Urine

Test Code : **8391808** - Laboratory Services

2

Method : **LC-MS/MS** - Performed : **Tue 08:00** - Reported : **2 weeks**

Sample Name	Container	Amount	Transfer
Urine, Frozen	Plastic Collection Cup	1mL	<0°C