



1	Troponin I-High Sensitive Test Code : 8349806 - Laboratory Services			
	Method : CLIA - Performed : Mon-Sat - Reported : 2 hours			
	Sample Name	Container	Amount	Transfer
	Serum	Separator Tube	4mL	2-8°C
	Plasma (Heparin)	Heparin-Yeşil Kapak(Jelli)	4mL	2-8°C

2	Cystatine-C Test Code : 8388301 - Laboratory Services			
	Method : NEPHELOMETRY - Performed : Mon 08:00 - Reported : 15 days			
	Sample Name	Container	Amount	Transfer
	Serum, Frozen	Separator Tube	1,5mL	<0°C

3	Cysteine-C, Urine Test Code : 8388303 - Laboratory Services			
	Method : NEPHELOMETRY - Performed : Mon 08:00 - Reported : 15 days			
	Sample Name	Container	Amount	Transfer
	Urine, Frozen	Plastic Collection Cup	2mL	<0°C

4	Complement C1q			
	Test Code : <b>8616011</b> - Laboratory Services			
	Method : <b>NEPHELOMETRY</b> - Performed : <b>Tue 08:00</b> - Reported : <b>1 weeks</b>			
	Sample Name	Container	Amount	Transfer
	Serum	Separator Tube	1mL	2-8°C
	Reference Ranges 			
Normal Values				
Gender	Lower Value	Upper Value	Unit	
	15.7	30.6	mg/dL	

5	Insulin Antibody			
	Test Code : <b>8633701</b> - Laboratory Services			
	Method : <b>IRMA</b> - Performed : <b>Mon 08:00</b> - Reported : <b>Tue 18:00</b>			
	Sample Name	Container	Amount	Transfer
Serum	Separator Tube	0,5mL	2-8°C	

6	Acetylcholine Receptor Antibody			
	Test Code : <b>8625610</b> - Laboratory Services			
	Method : <b>RIA</b> - Performed : <b>Mon 08:00</b> - Reported : <b>Tue 18:00</b>			
	Sample Name	Container	Amount	Transfer
Serum	Separator Tube	1mL	2-8°C	

7	<b>Plasma Renin Activity (PRA)</b> Test Code : <b>8424402</b> - Laboratory Services			
	Method : <b>RIA</b> - Performed : <b>Tue 08:00 Wed 08:00 Thu 08:00</b> - Reported : <b>Fri 18:00</b>			
	Sample Name	Container	Amount	Transfer
	Plasma (EDTA), Frozen	Lavender (EDTA)	1mL	<0°C

8	<b>Antidiuretic Hormone (ADH)</b> Test Code : <b>8458801</b> - Laboratory Services			
	Method : <b>RIA</b> - Performed : <b>Mon 08:00</b> - Reported : <b>Thu 18:00</b>			
	Sample Name	Container	Amount	Transfer
	Plasma (EDTA), Frozen	Lavender (EDTA)	3mL	<0°C

9	<b>Testosteron, Free</b> Test Code : <b>8440201</b> - Laboratory Services			
	Method : <b>RIA</b> - Performed : <b>Mon-Fri 08:00</b> - Reported : <b>Mon-Fri 18:00</b>			
	Sample Name	Container	Amount	Transfer
	Serum	Separator Tube	1mL	2-8°C

10

ALEX, Multiplex Allergy Test Panel

Test Code : **8019703** - Laboratory Services

Method : **EIA** - Performed : **Fri 08:00** - Reported : **Mon 18:00**

Sample Name	Container	Amount	Transfer
Serum	Separator Tube	0,5mL	2-8°C