

## Laboratory Test Updates

(effective on March 16, 2026)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8730102	Adenovirus Antigen, Stool		X								
2	C8730103	Adenovirus Antigen, Respiratory Specimen		X								
3	C8727001	Chlamydia trachomatis Antigen		X								
4	C8740001	Influenza Antigen (A and B)		X								
5	C8648001	QuantiFERON TB-Gold (Gamma Interferon Release Test) (IGRA)		X								
6	C8742001	Respiratory Syncytial Virus (RSV) Antigen		X								
7	C8742501	Rotavirus Antigen		X								
8	C8392908	SARS-CoV-2 Antigen (Rapid Test)		X								
9	C8788001	Streptococcus (Group A) Antigen, Rapid Test		X								
10	C8925714	Carbapenemase Gene Typing Test		X								
11	C8506001	Blood Smear Examination										X
12	C8430501	Somatomedin-C (IGF-1)					X					

1	<b>Adenovirus Antigen, Stool</b> Test Code : 8730102 - Laboratory Services			
	Methodology : IMMUNOCHROMATOGRAPHY - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
Stool	Plastic Collection Cup	1gr	2-8°C	

2	<b>Adenovirus Antigen, Respiratory Specimen</b> Test Code : 8730103 - Laboratory Services			
	Methodology : IMMUNOCHROMATOGRAPHY - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Nasopharyngeal Swab	Dacron/ Rayon Swab		2-8°C
	Nasopharyngeal Aspirate	Biologic Fluid Tube	2mL	2-8°C
	Transtracheal Aspirate	Biologic Fluid Tube	1mL	2-8°C
Bronchoalveolar Lavage	Biologic Fluid Tube	1mL	2-8°C	

3	<b>Chlamydia trachomatis Antigen</b> Test Code : 8727001 - Laboratory Services			
	Methodology : IMMUNOCHROMATOGRAPHY - Performed : Mon-Fri 08:00 Sat 08:00 - Reported : Tue-Sat 18:00 Mon 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Conjunctival Swab	Dacron/ Rayon Swab		2-8°C
	Urine	Plastic Collection Cup	20mL	2-8°C
	Cervical Swab	Dacron/ Rayon Swab		2-8°C
	Urethral Discharge Material	Dacron/ Rayon Swab		2-8°C
Prostatic fluid	Dacron/ Rayon Swab		2-8°C	

4	<b>Influenza Antigen (A and B)</b>		
	Test Code : 8740001 - Laboratory Services		
	Methodology : IMMUNOCHROMATOGRAPHY - Performed : Mon-Sat - Reported : 2 hours		
	Specimen Type	Collection Container/Tube	Specimen Volume
	Throat Swab	Dacron/ Rayon Swab	2-8°C
	Nasal Discharge	Dacron/ Rayon Swab	2-8°C
	Nasopharyngeal Aspirate	Dacron/ Rayon Swab	2-8°C
Nasopharyngeal Swab	Dacron/ Rayon Swab	2-8°C	
Sputum	Plastic Collection Cup	2-8°C	

5	<b>QuantiFERON TB-Gold (Gamma Interferon Release Test) (IGRA)</b>		
	Test Code : 8648001 - Laboratory Services		
	Methodology : EIA - Performed : Mon-Thu 08:00 Fri 08:00 - Reported : Wed-Sat 18:00 Mon 18:00		
Specimen Type	Collection Container/Tube	Specimen Volume	Transport
Contact to Laboratory	Call the Laboratory		2-8°C

6	<b>Respiratory Syncytial Virus (RSV) Antigen</b>		
	Test Code : 8742001 - Laboratory Services		
	Methodology : IMMUNOCHROMATOGRAPHY - Performed : Mon-Fri 12:00 Sat 12:00 - Reported : Tue-Sat 18:00 Mon 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
	Bronchoalveolar Lavage	Biologic Fluid Tube	1mL
	Transtracheal Aspirate	Biologic Fluid Tube	1mL
Nasopharyngeal Swab	Dacron/ Rayon Swab	2-8°C	
Nasopharyngeal Aspirate	Biologic Fluid Tube	2mL	

7	<b>Rotavirus Antigen</b>		
	Test Code : 8742501 - Laboratory Services		
	Methodology : IMMUNOCHROMATOGRAPHY - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00		
Specimen Type	Collection Container/Tube	Specimen Volume	Transport
Stool	Plastic Collection Cup	2gr	2-8°C

8	<b>SARS-CoV-2 Antigen (Rapid Test)</b>			
	Test Code : 8392908 - Laboratory Services			
	Methodology : IMMUNOCHROMATOGRAPHY - Performed : Mon-Sun - Reported : 1 hour			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
Nasopharyngeal Swab	Call the Laboratory		2-8°C	
Gargle and mouthwash	MyMagiCon Tube	20mL	18-24°C	

9	<b>Streptococcus (Group A) Antigen, Rapid Test</b>			
	Test Code : 8788001 - Laboratory Services			
	Methodology : IMMUNOCHROMATOGRAPHY - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00			
Specimen Type	Collection Container/Tube	Specimen Volume	Transport	
Throat Swab	Dacron/ Rayon Swab		18-24°C	

10	<b>Carbapenemase Gene Typing Test</b>			
	Test Code : 8925714 - Laboratory Services			
	Methodology : IMMUNOCHROMATOGRAPHY - Performed : Mon-Fri 08:00 - Reported : Mon-Fri 18:00			
Specimen Type	Collection Container/Tube	Specimen Volume	Transport	
Young Bacteria Strain	Plate		18-24°C	

12	<b>Somatomedin-C (IGF-1)</b>			
	Test Code : 8430501 - Laboratory Services			
	Methodology : CLIA - Performed : Tue 08:00 Thu 08:00 - Reported : Wed 18:00 Fri 18:00			
Specimen Type	Collection Container/Tube	Specimen Volume	Transport	
Serum, Frozen	Separator Tube	1mL	<0°C	