

Alzheimer panel, CSF

Test Code : **8217518** - Laboratory Services

1

Methodology : **ELISA** - Performed : **Tue 08:00** - Reported : **Fri 18:00**

Specimen Type	Container	Volume	Transport
CSF, Frozen	Polyethylene/Polypropylene Tube	1mL	<0°C

Additional Information

Transport: Frozen in polyethylene/polypropylene tube.

CSF, Frozen

Profile contains: Beta Amyloid (1-40), CSF, Beta Amyloid (1-42), CSF, Tau-protein, Phospho-Tau parameters.

Plasminogen Activity

Test Code : **8542103** - Laboratory Services

6

Methodology : **CHROMOGENIC** - Performed : **Tue 08:00** - Reported : **1 weeks**

Specimen Type	Container	Volume	Transport
Plasma (Citrate), Frozen	Blue (Sodium citrate)	1mL	<0°C

Additional Information

The plasma sample should be frozen by portioning it into at least two separate tubes.

Plasma (Citrate), Frozen

Unacceptable Conditions: Hemolyzed, icteric and lipemic specimens. Evaluating fibrinolysis, in combination with other components of the fibrinolytic system (fibrinogen, tissue plasminogen-activator-inhibitor, and d-dimers).