

Laboratory Test Updates

(effective on July 25, 2022)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	Update In Additional	New Tests	Internalized Tests	Reactivated Tests	Inactive
1	C8237401	Bicarbonate				X		X				
2	C8524002	ADAM TS 13 (von Willebrand Factor Splitting Protease)			X							
3	C8660306	Adenovirus Antibody, IgM			X							

1	Bicarbonate			
	Test Code : 8237401 - Laboratory Services			
	Method : PHOTOMETRY - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00			
	Sample Name	Container	Amount	Transfer
	Serum	Separator Tube	0,5mL	2-8°C
Additional Information ▼				
<p>Serum</p> <p>It is recommended to take the blood sample into a gel tube, centrifuge within 2 hours, and transfer it without opening the tube cap after centrifugation.</p>				

2

ADAM TS 13 (von Willebrand Factor Splitting Protease)

Test Code : **8524002** - Laboratory Services

Method : **ELISA** - Performed : **Mon 08:00** - Reported : **Wed 18:00**

Sample Name	Container	Amount	Transfer
Plasma (Citrate), Frozen	Blue (Sodium citrate)	4mL	<0°C

3

Adenovirus Antibody, IgM

Test Code : **8660306** - Laboratory Services

Method : **EIA** - Performed : **Fri 08:00** - Reported : **Mon 18:00**

Sample Name	Container	Amount	Transfer
Serum	Separator Tube	1mL	2-8°C

Reference Ranges

