

Laboratory Test Updates

(effective on June 30, 2025)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8252504	Copper (Cu), erythrocyte			X							
2	C8463002	Zinc (Zn), Erythrocyte			X							
3	C8373506	Magnesium (Mg), Erythrocyte			X							
4	C8425505	Selenium, RBCs			X							
5	C8666301	Anti-Epstein-Barr Virus (EBV)-EA Antibody			X							
6	C8660906	Borrelia burgdorferi Antibody, IgG (Lyme Disease)			X							
7	C8661802	Borrelia burgdorferi Antibody, IgM (Lyme Disease)			X							
8	C8674701	Parvovirus B19 Antibody, IgG			X							
9	C8674702	Parvovirus B19 Antibody, IgM			X							
10	C8673501	Mumps Antibody, IgG			X							
11	C8673502	Mumps Antibody, IgM			X							
12	C8676502	Rubeola (Measles) Antibody, IgG			X							
13	C8676501	Rubeola (Measles) Antibody, IgM			X							
14	C8678701	Varicella Zoster Virus (VZV) Antibody, IgG			X							
15	C8678702	Varicella Zoster Virus (VZV) Antibody, IgM			X							

1	Copper (Cu), erythrocyte		
	Test Code : 8252504 - Laboratory Services		
	Method : ICP-MS - Performed : Tue 08:00 Wed 08:00 Fri 08:00 - Reported : 1 week 18:00		
	Sample Name	Collection Container/Tube	Specimen Volume
	Whole Blood (EDTA)	Lavender (EDTA)	2mL
			Transfer
			2-8°C

2	Zinc (Zn), Erythrocyte		
	Test Code : 8463002 - Laboratory Services		
	Method : ICP-MS - Performed : Tue 08:00 Wed 08:00 Fri 08:00 - Reported : 1 week 18:00		
	Sample Name	Collection Container/Tube	Specimen Volume
	Whole Blood (EDTA)	Lavender (EDTA)	2mL
			Transfer
			2-8°C

3	Magnesium (Mg), Erythrocyte		
	Test Code : 8373506 - Laboratory Services		
	Method : ICP-MS - Performed : Tue 08:00 Wed 08:00 Fri 08:00 - Reported : 1 week 18:00		
	Sample Name	Collection Container/Tube	Specimen Volume
	Whole Blood (EDTA)	Lavender (EDTA)	2mL
			Transfer
			2-8°C

4	Selenium, RBCs		
	Test Code : 8425505 - Laboratory Services		
	Method : ICP-MS - Performed : Tue 08:00 Wed 08:00 Fri 08:00 - Reported : 1 week 18:00		
	Sample Name	Collection Container/Tube	Specimen Volume
	Whole Blood (EDTA)	Lavender (EDTA)	2mL
			Transfer
			2-8°C

5	Epstein-Barr Virus (EBV)-EA Antibody		
	Test Code : 8666301 - Laboratory Services		
	Method : CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00		
	Sample Name	Collection Container/Tube	Specimen Volume
	Serum	Separator Tube	1mL
			Transfer
			2-8°C

6	Borrelia burgdorferi Antibody, IgG (Lyme Disease) Test Code : 8660906 - Laboratory Services		
	Method : CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00		
	Sample Name	Collection Container/Tube	Specimen Volume
	Serum	Separator Tube	1mL

7	Borrelia burgdorferi Antibody, IgM (Lyme Disease) Test Code : 8661802 - Laboratory Services		
	Method : CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00		
	Sample Name	Collection Container/Tube	Specimen Volume
	Serum	Separator Tube	1mL

8	Parvovirus B19 Antibody, IgG Test Code : 8674701 - Laboratory Services		
	Method : CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00		
	Sample Name	Collection Container/Tube	Specimen Volume
	Serum	Separator Tube	1mL

9	Parvovirus B19 Antibody, IgM Test Code : 8674702 - Laboratory Services		
	Method : CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00		
	Sample Name	Collection Container/Tube	Specimen Volume
	Serum	Separator Tube	1mL

10	Mumps Antibody, IgG Test Code : 8673501 - Laboratory Services		
	Method : CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00		
	Sample Name	Collection Container/Tube	Specimen Volume
	Serum	Separator Tube	1mL

11	Mumps Antibody, IgM Test Code : 8673502 - Laboratory Services			
	Method : CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum	Separator Tube	1mL	2-8°C

12	Rubeola (Measles) Antibody, IgG Test Code : 8676502 - Laboratory Services			
	Method : CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum	Separator Tube	1mL	2-8°C

13	Rubeola (Measles) Antibody, IgM Test Code : 8676501 - Laboratory Services			
	Method : CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum	Separator Tube	1mL	2-8°C

14	Varicella Zoster Virus (VZV) Antibody, IgG Test Code : 8678701 - Laboratory Services			
	Method : CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum	Separator Tube	1mL	2-8°C

15	Varicella Zoster Virus (VZV) Antibody, IgM Test Code : 8678702 - Laboratory Services			
	Method : CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum	Separator Tube	1mL	2-8°C