



2	<b>Lupus Erythematosus (LE) cell</b>			
	Test Code : <b>8685002</b> - Laboratory Services			
	Method : <b>MICROSCOPY</b> - Performed : <b>Mon 08:00 Wed 08:00 Fri 08:00</b> - Reported : <b>Tue 18:00 Thu 18:00 Mon 18:00</b>			
	Sample Name	Container	Amount	Transfer
Serum	Separator Tube	0,5mL	2-8°C	

3	<b>Echinococcus Antibody, IgG</b>			
	Test Code : <b>8668203</b> - Laboratory Services			
	Method : <b>EIA</b> - Performed : <b>Tue 08:00</b> - Reported : <b>Wed 18:00</b>			
	Sample Name	Container	Amount	Transfer
Serum	Separator Tube	0,3mL	2-8°C	

4	<b>Echinococcus Antibody, IgM</b>			
	Test Code : <b>8668204</b> - Laboratory Services			
	Method : <b>EIA</b> - Performed : <b>Tue 08:00</b> - Reported : <b>Wed 18:00</b>			
	Sample Name	Container	Amount	Transfer
Serum	Separator Tube	0,3mL	2-8°C	