

# Laboratory Test Updates

(effective on July 04, 2022)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	Update In Additional	New Tests	Internalized Tests	Reactivated Tests	Inactive
1	C8780236	Monkeypox Virus DNA PCR							X			
2	C8600418	Specific IgE- Vanilla (f234)								X		
3	C8210801	Aluminum				X						
4	C8301813	Antimony	X			X						
5	C8252501	Copper (Cu)				X						
6	C8463001	Zinc (Zn)				X						
7	C8301801	Cobalt	X			X						
8	C8378501	Manganese	X			X						
9	C8301826	Molybdenum	X			X						
10	C8388501	Nickel				X						
11	C8425501	Selenium				X						
12	C8412007	Vanadium				X						
13	C8382501	Mercury				X						
14	C8210803	Tin				X						
15	C8301827	Molybdenum, Whole Blood	X			X						

1	<b>Monkeypox Virus DNA PCR</b>			
	Test Code : <b>8780236</b> - Laboratory Services			
	Method : <b>PCR</b> - Performed : <b>Mon-Fri</b> - Reported : <b>1 day</b>			
	<b>Sample Name</b>	<b>Container</b>	<b>Amount</b>	<b>Transfer</b>
	Lesion material	Nucliswab		2-8°C

2	<b>Specific IgE- Vanilla (f234)</b>			
	Test Code : <b>8600418</b> - Laboratory Services			
	Method : <b>CLIA</b> - Performed : <b>Mon-Thu 08:00 Fri 08:00</b> - Reported : <b>Tue-Fri 18:00 Mon 18:00</b>			
	<b>Sample Name</b>	<b>Container</b>	<b>Amount</b>	<b>Transfer</b>
	Serum	Separator Tube	0,5mL	2-8°C

3	<b>Aluminum</b>			
	Test Code : <b>8210801</b> - Laboratory Services			
	Method : <b>ICP-MS</b> - Performed : <b>Tue 08:00 Wed 08:00 Thu 08:00</b> - Reported : <b>1 week 18:00</b>			
	<b>Sample Name</b>	<b>Container</b>	<b>Amount</b>	<b>Transfer</b>
	Plasma (Heparin)	Sodium Heparin-Royal Blue	3mL	2-8°C

4	<b>Antimony</b>			
	Test Code : <b>8301813</b> - Laboratory Services			
	Method : <b>ICP-MS</b> - Performed : <b>Tue 08:00 Wed 08:00 Thu 08:00</b> - Reported : <b>1 week 18:00</b>			
	<b>Sample Name</b>	<b>Container</b>	<b>Amount</b>	<b>Transfer</b>
	Plasma (Heparin)	Sodium Heparin-Royal Blue	2mL	2-8°C

5	Copper (Cu)			
	Test Code : <b>8252501</b> - Laboratory Services			
	Method : <b>AAS</b> - Performed : <b>Tue 08:00 Thu 08:00</b> - Reported : <b>Wed 18:00 Fri 18:00</b>			
	Sample Name	Container	Amount	Transfer
Plasma (Heparin)	Sodium Heparin-Royal Blue	2mL	2-8°C	

6	Zinc (Zn)			
	Test Code : <b>8463001</b> - Laboratory Services			
	Method : <b>AAS</b> - Performed : <b>Mon 08:00 Wed 08:00 Fri 08:00</b> - Reported : <b>Tue 18:00 Thu 18:00 Mon 18:00</b>			
	Sample Name	Container	Amount	Transfer
Plasma (Heparin)	Sodium Heparin-Royal Blue	2mL	2-8°C	

7	Cobalt			
	Test Code : <b>8301801</b> - Laboratory Services			
	Method : <b>ICP-MS</b> - Performed : <b>Tue 08:00 Wed 08:00 Thu 08:00</b> - Reported : <b>1 week 18:00</b>			
	Sample Name	Container	Amount	Transfer
Plasma (Heparin)	Sodium Heparin-Royal Blue	2mL	2-8°C	

8	Manganese			
	Test Code : <b>8378501</b> - Laboratory Services			
	Method : <b>ICP-MS</b> - Performed : <b>Tue 08:00 Wed 08:00 Thu 08:00</b> - Reported : <b>1 weeks 18:00</b>			
	Sample Name	Container	Amount	Transfer
Plasma (Heparin)	Sodium Heparin-Royal Blue	2mL	2-8°C	

9	<b>Molybdenum</b> Test Code : <b>8301826</b> - Laboratory Services			
	Method : <b>ICP-MS</b> - Performed : <b>Tue 08:00 Wed 08:00 Thu 08:00</b> - Reported : <b>1 weeks 18:00</b>			
	<b>Sample Name</b>	<b>Container</b>	<b>Amount</b>	<b>Transfer</b>
	Plasma (Heparin)	Sodium Heparin-Royal Blue	2mL	2-8°C

10	<b>Nickel</b> Test Code : <b>8388501</b> - Laboratory Services			
	Method : <b>ICP-MS</b> - Performed : <b>Tue 08:00 Wed 08:00 Thu 08:00</b> - Reported : <b>1 weeks 18:00</b>			
	<b>Sample Name</b>	<b>Container</b>	<b>Amount</b>	<b>Transfer</b>
	Plasma (Heparin)	Sodium Heparin-Royal Blue	3mL	2-8°C

11	<b>Selenium</b> Test Code : <b>8425501</b> - Laboratory Services			
	Method : <b>ICP-MS</b> - Performed : <b>Tue 08:00 Wed 08:00 Thu 08:00</b> - Reported : <b>1 week 18:00</b>			
	<b>Sample Name</b>	<b>Container</b>	<b>Amount</b>	<b>Transfer</b>
	Plasma (Heparin)	Sodium Heparin-Royal Blue	2mL	2-8°C

12	<b>Vanadium</b> Test Code : <b>8412007</b> - Laboratory Services			
	Method : <b>ICP-MS</b> - Performed : <b>Tue 08:00 Wed 08:00 Thu 08:00</b> - Reported : <b>1 week 18:00</b>			
	<b>Sample Name</b>	<b>Container</b>	<b>Amount</b>	<b>Transfer</b>
	Plasma (Heparin)	Sodium Heparin-Royal Blue	2mL	2-8°C

13	<b>Mercury</b> Test Code : <b>8382501</b> - Laboratory Services			
	Method : <b>ICP-MS</b> - Performed : <b>Tue 08:00 Wed 08:00 Thu 08:00</b> - Reported : <b>1 week 18:00</b>			
	<b>Sample Name</b>	<b>Container</b>	<b>Amount</b>	<b>Transfer</b>
	Whole Blood (Heparin)	Sodium Heparin-Royal Blue	2mL	2-8°C

14	Tin			
	Test Code : <b>8210803</b> - Laboratory Services			
	Method : <b>ICP-MS</b> - Performed : <b>Tue 08:00 Wed 08:00 Thu 08:00</b> - Reported : <b>1 weeks 18:00</b>			
	<b>Sample Name</b>	<b>Container</b>	<b>Amount</b>	<b>Transfer</b>
	Plasma (Heparin)	Sodium Heparin-Royal Blue	3mL	2-8°C

15	Molybdenum, Whole Blood			
	Test Code : <b>8301827</b> - Laboratory Services			
	Method : <b>ICP-MS</b> - Performed : <b>Tue 08:00 Wed 08:00 Thu 08:00</b> - Reported : <b>1 weeks 18:00</b>			
	<b>Sample Name</b>	<b>Container</b>	<b>Amount</b>	<b>Transfer</b>
	Whole Blood (Heparin)	Sodium Heparin-Royal Blue	5mL	2-8°C